


T MAZE



Note: Drawing used for illustration purposes.

Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

T MAZE SPECIFICATIONS

MAZE DIMENSIONS

Mouse	(A) 24"	(B) 3"	(C) 10.5"	(D) 15"	(E) 5"	(F) 3.5"	(H) .75" Diam. x .5" Deep
Rat	(A) 36"	(B) 4"	(C) 16"	(D) 22"	(E) 10"	(F) 7"	(H) 1" Diam. x .75" Deep
Material Composition	1/4" Haircell ABS Plastic						
Color Options	Black, White, Beige						

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
T Maze for Mice	Black 7001-0404	T Maze for Rats	Black 7001-0413
	White 7001-0401		White 7001-0410
	Beige 7001-0407		Beige 7001-0416

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.