


Y MAZE



Note: Drawing used for illustration purposes.

Features & Benefits

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

Y MAZE SPECIFICATIONS

MAZE DIMENSIONS

Mouse	(A) 3"	(B) 15"	(C) 15"	(D) 5"	(E) 3.5"	(H) .75" Diam. X .5" Deep
Rat	(A) 4"	(B) 22"	(C) 22"	(D) 10"	(E) 7"	(H) 1" Diam. X .75" Deep
Material Composition	1/4" Haircell ABS Plastic					
Color Options	Black, White, Beige					

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Y Maze for Mice	Black 7001-0422	Y Maze for Rats	Black 7001-0431
	White 7001-0419		White 7001-0428
	Beige 7001-0425		Beige 7001-0434

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.