


BARNES MAZE



Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One target box and 19 false boxes
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

BARNES MAZE SPECIFICATIONS

Maze Dimensions & No. of Holes	Rat: 48" diameter with 20 - 4" holes Mouse: 36" diameter with 20 - 2" holes
Stand	36" (H)
Start Tube	Rat: 12"(H) , 8"(D) Mouse: 8" (H), 6"(D)
Weight	48" (Rat Model): 48 lbs. 36" (Mouse Model): 25 lbs.
Material Composition	Top: White, beige or black plastic Target & False Target Boxes: Black plastic

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Barnes Maze for Mice	Black 7001-0324	Barnes Maze for Rats	Black 7001-0325
	White 7001-0235		White 7001-0326
	Beige 7001-0340		Beige 7001-0341

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.