


MULTIPLE UNIT OPEN FIELD ENCLOSURE



Features & Benefits

- Made of high density, non-porous plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

MULTIPLE UNIT OPEN FIELD ENCLOSURE SPECIFICATIONS

Base	3 ft. 6" x 3 ft. 6"
Outer Wall Length	3 ft. 4" x 3 ft. 4"
Wall Height (Inner & Outer)	1 ft. 3"
Inside Chamber Footprint	1 ft. 7 1/2" x 1 ft. 7 1/2"
Color Options	Black, White, Beige

Black 7001-0353

White 7001-0354

Beige 7001-0067

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.