


RADIAL ARM MAZE



Note: Drawing used for illustration purposes.

Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

RADIAL ARM MAZE SPECIFICATIONS

	MAZE DIMENSIONS					
Mouse	(A) 26.45"	(B) 9"	(C) 8.45"	(D) 3"	(E) 6"	(H) .75" Diam. x .5" Deep
Rat	(A) 47.86"	(B) 18.5"	(C) 10.86"	(D) 4"	(E) 10"	(H) 1" Diam. x .75" Deep
Material Composition	1/4" ABS Plastic					
Color Options	Black, White, Beige					

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Radial Arm Maze for Mice	Black 7001-0440	Radial Arm Maze for Rats	Black 7001-0449
	White 7001-0437		White 7001-0446
	Beige 7001-0443		Beige 7001-0452

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.