


RADIAL ARM MAZE



Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI video tracking systems that feature:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

Note: Drawing used for illustration purposes.

RADIAL ARM MAZE SPECIFICATIONS

	MAZE DIMENSIONS				
Mouse	(A) 24"	(B) 9"	(C) 5.25"	(D) 2"	(E) 6"
Rat	(A) 48"	(B) 18.5"	(C) 11.5"	(D) 4"	(E) 8"
Material Composition	1/4" ABS Plastic				
Color Options	Black, White, Beige				

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Radial Arm Maze for Mice	Black 7001-0186	Radial Arm Maze for Rats	Black 7001-0028
	White 7001-0343		White 7001-0328
	Beige 7001-0342		Beige 7001-0344

NOTE: Beige is SDI's standard Maze color unless specified by the researcher. The SMART Video Tracking System Software is sensitive enough to track a subject with little or no contrast.