



WATER MAZE SET UP COMPONENTS

Features

- Tank available in blue or white
- Platform base and post clear acrylic. Attached top available in clear, white or black acrylic
- Shipping included in the price of the tank
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports



Note: Drawings and images are not shown proportionally for illustration purposes.

WATER MAZE SET UP COMPONENTS SPECIFICATIONS

	TANK DIMENSIONS			PLATFORM DIMENSIONS		
Mouse	(A) 48"	(B) 30"	210 GAL, 49 LBS (empty)	Square (C) 12" x 12"	(D) 12.76"	(E) 4" x 4"
Rat	(A) 72"	(B) 30"	500 GAL, 73 LBS (empty)	Round (C) 12" x 12"	(D) 12.76"	(F) 4" diameter
Platform Material Composition	Acrylic					
Tank Material Composition	Polyethylene					
Tank Color Options	Blue, White					

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Water Tank for Mice (4')	Blue 7000-0720	Water Tank for Rat (6')	Blue 7000-0718
	White 7000-0723		White 7000-0719
Water Maze Platforms- Square Tops	Clear 7500-0271	Water Maze Platforms- Round Tops	Clear 7500-0274
	White 7500-0272		White 7500-0275
	Black 7500-0273		Black 7500-0276
Drain Kit	7500-0217		