

## WHITE NOISE GENERATOR

*Programmable and manual scrambled shock administration in one solid state unit*



### PRODUCT OVERVIEW


SDI's White Noise Generator generates adjustable amplitude, constant white noise to eliminate extraneous background sounds. The White Noise Generator is a compact, solid state unit that is small and lightweight enough to mount on a laboratory wall. It can also be installed directly into an animal testing enclosure. The White Noise Generator blocks out other stimuli that may cause distraction during behavior testing.

### WHITE NOISE GENERATOR COMPONENTS

- White Noise Generator Unit
- Power Supply

### Features & Benefits

- Compact, solid state unit
- Lightweight for easy mounting on a wall
- Can be directly installed in a test enclosure



## WHITE NOISE SPECIFICATIONS

Dimensions	3 ¼" (W) x 4 ¾" (D) x 2 ¼" (H)
Weight	1 lb.
Material Composition	Gray plastic
Standard Cable Length	6 ft.

## FOR MORE INFORMATION

To learn more about SDI behavioral testing systems, please visit [www.sandiegoinstruments.com](http://www.sandiegoinstruments.com). If you have any questions or would like to request a quote please call (858) 530-2600 or email us at [sales@sandiegoinstruments.com](mailto:sales@sandiegoinstruments.com).

## SDI TEST ACCESSORIES

- AIRSTIM™
- AIRSTIM™ with software
- Portable Air Puff Unit
- Portable Air Puff Unit with software
- Animal Shocker
- Animal Stimulator



San Diego Instruments, Inc.  
6295 Ferris Square, Suite A  
San Diego, CA 92121  
Ph: 858-530-2600  
Fax: 858-530-2646  
[www.sandiegoinstruments.com](http://www.sandiegoinstruments.com)

©2010 San Diego Instruments. All rights reserved. SDI and the SDI logo are trademarks of San Diego Instruments, Inc. All other trademarks mentioned herein are property of their respective owners. Specifications are subject to change without notice. The equipment described herein is designed for research and educational purposes and is not intended for the diagnosis, alleviation, treatment, monitoring or prevention of disease, injury or handicap.  
0910-whitenoise