


ZERO MAZE



Note: Drawing used for illustration purposes.

Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI ANY-Maze™ video tracking system that features:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

ZERO MAZE SPECIFICATIONS

	MAZE DIMENSIONS			3 Adjustable Legs (not shown)
	Inside Wall	Outside Wall	(C)	
Mouse	(A) 20" DIA. / 62.83" CIRCUM.	(B) 24" DIA. / 75.40" CIRCUM.	(C) 6"	24"
Rat	(A) 40" DIA. / 125.66" CIRCUM.	(B) 48" DIA. / 150.8" CIRCUM.	(C) 8"	30"
Material Composition	1/4" ABS Plastic			
Color Options	Black, White, Beige			

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Zero Maze for Mice	Black 2325-0230	Zero Maze for Rats	Black 2325-0227
	White 2325-0229		White 2325-0226
	Beige 2325-0231		Beige 2325-0228

NOTE: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-Maze™ Video Tracking System.